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Introduction to immunogenicity assessment for  
protein therapeutics

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(Introduction to immunogenicity assessment for protein therapeutics)

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Some factors which influence the unwanted immunogenicity of therapeutic proteins.	
Factors determining the immunogenicity of therapeutic proteins	
Product-related (Intrinsic)	
<ul style="list-style-type: none"> <li>• molecular structure – primary amino acid sequence or variants, novel epitopes;</li> <li>• aggregates, degradation products, oxidized or deamidated forms;</li> <li>• host cell DNA/proteins;</li> <li>• formulation;</li> <li>• primary packaging.</li> </ul>	
Product-related (Extrinsic)	
<ul style="list-style-type: none"> <li>• dose, route, frequency of administration, episodic/continuous treatment, duration of treatment, previous exposure;</li> <li>• cellular or soluble target;</li> <li>• biological properties of the therapeutic;</li> <li>• endogenous counterpart, redundant or non-redundant.</li> </ul>	
Host-related	
<ul style="list-style-type: none"> <li>• genetic profile;</li> <li>• immune status, disease state, medication.</li> </ul>	

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Commonly used screening assays.		
Type of Assay	Advantages	Disadvantages
Direct/Indirect ELISA	Easy to use and automate High through-put High therapeutic tolerance Inexpensive Generic reagents and instrument	May bind non-specifically High background May fail to detect low-affinity antibodies Requires species specific secondary reagent
Bridging ELISA	Easy to use and automate High through-put Low background, High therapeutic tolerance in solution phase High specificity (dual-arm binding) Generic reagents and instrument	Antigen labelling required May fail to detect low-affinity antibodies Highly susceptible to interference by therapeutic, serum components e.g., anti-human Ig molecules, multivalent targets May not detect IgG4
Electrochemiluminescence	High through-put, large dynamic range Minimally affected by matrix High tolerance to therapeutic in solution phase Detection signal consistent during life of TAG conjugate	May require two antigen conjugates Antigen labelling required Susceptible to interference by therapeutic, serum components e.g., anti-human Ig molecules, multivalent targets May not detect IgG4
Radioimmunoprecipitation assay	Moderate through-put High sensitivity Can be specific Inexpensive	Vendor-specific equipment & reagents Can be isotype specific May not detect low-affinity antibodies. Requires radiolabelled antigen. Decay of radio-label may affect antigen stability
Surface plasmon resonance	Automated Determines specificity, isotype, relative binding affinity Enables detection of both 'low-affinity' and high affinity antibodies. Detection reagent not required	Antigen immobilization may alter therapeutic. Regeneration step may degrade antigen. Sensitivity often less than binding assay. Expensive vendor-specific equipment & reagents

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1. World Health Organization. Guidelines on evaluation of similar biotherapeutic products (SBPS).
2. European Medicines Agency. Guideline on immunogenicity assessment of biotechnology-derived therapeutic proteins.
3. FDA Guidance for Industry. Assay development and validation for immunogenicity testing of therapeutic proteins.
4. FDA Guidance for Industry. Immunogenicity Assessment for Therapeutic Protein Products.
5. M. Wadhwa, I. Knezevic, H. N. Kang, R. Thorpe. Immunogenicity assessment of biotherapeutic products: An overview of assays and their utility. *Biologicals*. 2015, 43: 298-306.
6. C. Pineda, G. C. Hernández, I. A. Jacobs, D. F. Alvarez, C. Carini. *BioDrugs*. Assessing the Immunogenicity of Biopharmaceuticals. *BioDrugs*. 2016, 30: 195-206.
7. A. Mikulskis, D. Yeung, M. Subramanyam, L. Amaravadi. Solution ELISA as a platform of choice for development of robust, drug tolerant immunogenicity assays in support of drug development. *J Immunol Meth*. 2010, 365 (1-2). 38-49.
8. J.M. Carrascosa. Immunogenicity in Biologic Therapy: Implications for Dermatology. *Actas Dermosifiliogr*. 2013, 104 (6): 471-479.
9. M. Spengler, M. Adlera, C. M. Niemeyer. Highly sensitive ligand-binding assays in pre-clinical and clinical applications: immuno-PCR and other emerging techniques. *Analyst*. 2015, 140 (18):6175-6194.
10. S. Mariani, M. Minunni. Surface plasmon resonance applications in clinical analysis. *Anal Bioanal Chem*. 2014, 406 (9-10): 2303-2323.
11. J. A. Lofgren, S. Dhandapani, J. J. Pennucci, C. M. Abbott, D. T. Mytych, A. Kaliyaperumal, S. J. Swanson, M. C. Mullenix. Comparing ELISA and Surface Plasmon Resonance for Assessing Clinical Immunogenicity of Panitumumab. *J Immunol*. 2007, 178 (11): 7467-7472.
12. J. C. Brose, P. J. Couble, M. R. Deehan, R. J. Nelson, W. G. Ferlin, S. Lory. Evaluation of Multiple Immunoassay Technology Platforms to Select the Anti-Drug Antibody Assay Exhibiting the Most Appropriate Drug and Target Tolerance. *J Immunol Res*. 2016, 1-15.
13. P. A. van Schouwenburg, T. Rispens, G. J. Wolbink. Immunogenicity of anti-TNF biologic therapies for rheumatoid arthritis. *Nat Rev Rheumatol*. 2016, 16 (2): 201-211.
14. C. Núñez, M. C. Cé nit, R. Alvarez-Lafuente, J. R í o, M. Fern á ndez-Arquero, R. Arroyo, X. Montalb á n, O. Fern á ndez, B. Oliver-Martos, L. Leyva, M. Comabella, E. Urcelay. HLA alleles as biomarkers of high-titre neutralising antibodies to interferon- $\beta$  therapy in multiple sclerosis. *J Med Genet*. 2014, 51 (6): 395-400.

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## Expert Analytical Services for Biologics

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#### LC-MS Based Analysis for Biologics Development

- mAb, ADC, Peptides, Glycoproteins, Recombinant Proteins, Vaccine
- ☑ Full Structural Characterization
- ☑ Determination of Impurities
- ☑ Production Optimization
- ☑ Comparability Studies of Biosimilars
- ☑ QC Release Testing
- ☑ Stability Testing

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#### GLP Bioanalytical Services in Accordance with ICH Q6B

- ☑ Physicochemical Characterization
  - Molecular Weight and Size
  - Amino Acid Sequencing
  - Terminal Sequencing
  - Peptide Mapping
  - Disulfide Bridges
  - Carbohydrate Structure
- ☑ Process- and Product-related Impurities
  - Truncated Forms
  - Post-translational Modifications
  - Host Cell Proteins
  - Antibiotics

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#### Bioanalysis & Customized Method Development

- ☑ ELISA Assay
- ☑ BA/BE/PK Studies
- ☑ Biomarker Discovery and Verification
- ☑ Immunogenicity