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

Table 1. Overview of methods for HCP detection

	Method	Comment	Limit of detection	Type of host cell	
Immuno-specific	ELISA		1.0 ng/mL	<i>E. coli</i> , SP 2/0	
	Threshold		0.5 ng/mL	<i>E. coli</i>	
	Western blot		Colorimetric (e.g. BCIP/NBT ^a)	0.2–0.5 ng/protein	<i>E. coli</i>
			Chemiluminescent (e.g. TMB ^b)	0.1–0.3 ng/protein	SP 2/0
	Slot blot		6.0 ng/mL	<i>E. coli</i>	
	Automated PLA ^c		0.025 ng/mL	CHO	
	Biosensors		10 ⁻¹⁸ M	<i>E. coli</i>	
Non-specific	1D and 2D electrophoresis	Silver stain	2.0–13.0 ng/protein ^d	<i>E. coli</i> , SP 2/0	
		Coomassie	8.0–52.0 ng/protein ^d	<i>E. coli</i>	
		SYPRO Ruby	0.3–1.0 ng/protein ^d	<i>E. coli</i> , CHO	
		CyDyes	0.03–0.08 ng/protein ^d	<i>E. coli</i> , CHO	
		Autoradiography	0.1–2 pg/protein	CHO	
	MS		1–50 fmol	<i>E. coli</i> , CHO	
	2D chromatography	MS	50 ppm for complex samples	<i>E. coli</i>	

a) BCIP = 5-bromo-4-chloro-3-indolyl-phosphate, NBT = nitro blue tetrazolium. b) TMB = 3,3',5,5'-tetramethylbenzidine. c) PLA = proximity ligation assay. d) Depending on protein and staining protocol.

表1. HCP检测方法分類^[6](Biotechnol. J. 2013, 8: 655-70)。

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1. USP Monographs: Insulin Human.
2. D. G. Bracewell, R. Francis, C. M. Smales. The Future of Host Cell Protein (HCP) Identification During Process Development and Manufacturing Linked to a Risk-Based Management for Their Control. *Biotechnol. bioeng.* 2015, 112 (9): 1727-37.
3. X. Wang, A. K. Hunter, N. M. Mozier. Host Cell Proteins in Biologics Development: Identification, Quantitation and Risk Assessment. *Biotechnol. bioeng.* 2009, 103 (3): 446-58.
4. K. Hoffman. Strategies for Host Cell Protein Analysis. *Biopharm.* 2000, 13 (6): 38-45.
5. USP Chapter <1132>: Residual Host Cell Protein Measurement in Biopharmaceuticals
6. A. L. Tscheliessnig, J. Konrath, R. Bates, et al. Host cell protein analysis in therapeutic protein bioprocessing-methods and applications. *Biotechnol. J.* 2013, 8: 655-70.
7. J. Zhu-Shimoni, C. Yu, J. Nishihara, et al. Host Cell Protein Testing by ELISAs and the Use of Orthogonal Methods. *Biotechnol. bioeng.* 2014. 111 (12): 2367-79.
8. X. Lu, C. F. Xu, S. Liu, et al. Identification and Quantification of Heat-Shock Protein 70: A Major Host-Cell Protein Contaminant from HEK Host Cells. *Bioprocess Intl.* 2015.
9. N. Liu, M. axim Brevnov, M. Furtado, et al. Host cellular protein quantification. 2012. *Bioprocess Intl.* 10 (2): 44-50.
10. N. Dagouassat, J. F. Haeuw, V. Robillard, et al. Development of a quantitative assay for residual host cell proteins in a recombinant subunit vaccine against human respiratory syncytial virus. *J. Immunol. Methods.* 2001. 251 (1): 151-59.
11. K. Teeparuksapun, M. Hedström, P. Kanatharana, et al. Capacitive immunosensor for the detection of host cell proteins. *J. Biotechnol.* 2010, 157, 207–13.
12. T. Rabilloud, M. Chevallet, S. Luche, et al. Two-dimensional gel electrophoresis in proteomics: Past, present and future. *J. Proteomics* 2010, 73, 2064–77.
13. S. H. Stansfield, E. E. Allen, D. M. Dinnis, et al. Dynamic analysis of GS-NS0 cells producing a recombinant monoclonal antibody during fed-batch culture. *Biotechnol. Bioeng.* 2007, 97, 410–24.
14. J. Peng, J. E. Elias, C. C. Thoreen, et al. Evaluation of multidimensional chromatography coupled with tandem mass spectrometry (LC/LC-MS/MS) for large-scale protein analysis: the yeast proteome. *J. Proteome Res.* 2003, 2(1): 43-50.

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Expert Analytical Services for Biologics

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LC-MS Based Analysis for Biologics Development

- mAb, ADC, Peptides, Glycoproteins, Recombinant Proteins, Vaccine
- ☑ Full Structural Characterization
- ☑ Determination of Impurities
- ☑ Production Optimization
- ☑ Comparability Studies of Biosimilars
- ☑ QC Release Testing
- ☑ Stability Testing

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GLP Bioanalytical Services in Accordance with ICH Q6B

- ☑ Physicochemical Characterization
 - Molecular Weight and Size
 - Amino Acid Sequencing
 - Terminal Sequencing
 - Peptide Mapping
 - Disulfide Bridges
 - Carbohydrate Structure
- ☑ Process- and Product-related Impurities
 - Truncated Forms
 - Post-translational Modifications
 - Host Cell Proteins
 - Antibiotics

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Bioanalysis & Customized Method Development

- ☑ ELISA Assay
- ☑ BA/BE/PK Studies
- ☑ Biomarker Discovery and Verification
- ☑ Immunogenicity